


Switch Stream

Quick Start Guide

Switch Stream


Quick Start Guide

1. Power button*
2. USBC charging port
3. Speaker



Input Side


4. Four assistive inputs
5. Ring Light



Output Side

for 2 way monocord only

6. Four outputs for the two way mono cord



*To reset the hub, hold the Power button in for 15-20 seconds

In The Box

Hub	x1
USBC Power Cable	x1
Two way monocord	x1
MISO Battery Pack	x1
Quick Start Guide	x1



Need Help?

Tutorial videos www.youtube.com/switchstream

Email info@getswitchstream.com

Call 267-799-6135

Set up

- 1. Power on:** using the included USB-C cable, connect the Hub to a power source like your laptop or a wall outlet.
 - For portability, use the included MISO battery pack as your power source - plug the USBC cord into the Hub's USBC port.
- 2.** Go to your Phone's Bluetooth settings and select **Switch Stream Hub**. Accept the pairing request.
- 3.** You should have received an email from "TestFlight" inviting you to download the Switch Stream app. (If you did not receive the email please check your spam folder. If you still do not see the email, please email info@getswitchstream.com) You can also use the link on the bottom of the Hub to download the app.
- 4.** Follow the instructions in the Switch Stream App and tap "Continue Setup" to finish setup.

Using switches with your Hub

To use switches with the Hub, plug your switches into the input side. Configure the switches in the app.

Using switch adapted devices with your Hub

To use switch adapted devices with the Hub, **plug your switches into the input side** and **plug the two way monocord into the output side**. Then plug your adapted device into the two way mono cord. Configure the switch input and output in the app.

Communicator / Audio

To record a message:

1. On the home screen of the app, tap any Switch Input
2. Tap the red circle to start record your message, then tap it again to stop recording
3. Your new recording shows up at the bottom of the list
4. To name the recording, tap the recording and add a name where it says "New Name" then tap submit
5. To set the recording to your switch check the box next to your recording and tap "Send Recording to Hub". You can select up to 4 recordings to play on one switch.
6. Your switch is now set up to play your recording(s). The recording will play through the Hub's speaker.

Keyboard / Switch Control

To use Switch Control

- Plug in the Hub
- Make sure your phone's Bluetooth is **on**.
- In your phone's Bluetooth settings find '**Switch Steam Hub**' and accept the pairing request.
- Open the Switch Stream app and tap Continue setup.
- Tap 'Connect to Smart Hub' which will scan for the Hub, tap the hub that pops up.
- You are now ready to set your switches up.
- Go to the keyboard screen and tap a key or choose switch control.
- If you are using an external iPad/device for switch control, after setting up your switches in the app:
 - Minimize the app (tap the home button)
 - Go to your phone's Bluetooth settings and **disconnect** from 'Switch Stream Hub' or forget device.
 - In your external device's Bluetooth settings, find 'Switch Stream Hub' and **accept the pairing request**.
 - You are now ready to use Switch Control on that device
- When using Switch Control, select 'external switch' and follow on screen instructions

Troubleshooting steps if Switch Control / Keys are not working:

Go to **Settings > Accessibility > Keyboards > turn Full Keyboard Access on**

Switch Modes - for anything plugged into Output

To use a switch adapted toy / device

- Using the included two-way monocord, plug a switch adapted device into any hub output 1-4.
- In the app on the plugged in screen, use the Output spinner to select the output you plugged into
- Tap 'Add Task' for any of the following switch modes
 - **Direct Press:** hold the switch down to activate
 - **Latch on/off (1 switch):** activation 1 = On, activation 2 = Off
 - **Latch on (2 switch):** activation Switch 1 = On
 - **Latch off (2 switch, set a second input in app, select the corresponding output):** activation Switch 2 = Off
 - **Timed Latch:** 1 activation turns device on for selected time
- Your switch should now activate the switch adapted device in the corresponding switch mode